VBM

DE

NEU: Video Animation
Larynx-Tubus inklusive
Cuffdruckmessung
www.vbm-medical.de/cuffdruck

EN

NEW: Video Animation
Laryngeal Tube including
Cuff Pressure Measurement
www.vbm-medical.de/cuffpressure

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EN

Manufacturers Note Cuff Pressure Measurement

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EN

Dear user,

the Laryngeal Tube is an established supraglottic airway device that has been implemented into various international practice guidelines. The Laryngeal Tube distinguishes itself through an effective, fast and straightforward handling to secure the airway – especially in case of emergency. It is the only supraglottic airway device with two high volume low pressure cuffs to ensure a reliable seal whilst at the same time being atraumatic to the mucosa. However, this is only the case if the correct cuff pressure is used.

VBM would like to stress the importance of cuff pressure measuring as it is described in the VBM instruction manual:

Excessive cuff pressure can lead to swelling of the tongue.
Adjust the cuff pressure to maximum 60 cmH₂O via the VBM cuff pressure gauge.

Particularly in the prehospital setting, but also in routine use, it appears monitoring cuff pressure is being neglected for various reasons. This can lead to tongue swelling (Bernhard et al. 2014, Schalk et al. 2014) which can only be avoided by monitoring and adjusting cuff pressure (Schalk et al. 2015, Kriege et al. 2017). For this purpose VBM offers different cuff pressure gauges which reduce the risk of pressure necrosis and mucosal ischemia.

Moreover, assessing, controlling and monitoring the correct position and adequate ventilation of the Laryngeal Tube in the patient, should be integrated into the training of rescue personnel (cf. Bernhard et al. 2014).

If you have any questions regarding “cuff pressure”, please do not hesitate to contact us under cuffpressure@vbm-medical.de.

Order information

Cuff Manometer **Universal**

Analog Cuff Pressure Gauge, complete with connecting tube (100 cm).

Scala (Ø 68 mm) with two green ranges

- for tracheal tubes (22-32 cmH₂O)
- for Laryngeal Tubes and Laryngeal Masks (32-60 cmH₂O)

REF 54-07-000

References

Bernhard, M. et al.: Prehospital airway management using the laryngeal tube, *Der Anaesthesist* 63.7 (2014) 589-596.

Schalk, R. et al.: Complications associated with the prehospital use of laryngeal tubes, *Resuscitation* 85 (2014) 1629-1632.

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Kriege, M. et al.: Evaluation of the optimal cuff volume and cuff pressure of the revised laryngeal tube "LTS-D" in surgical patients, *BMC Anesthesiology* (2017) Vol. 17, p1-7.7p.

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VBM **Medizintechnik GmbH**